

Chest freezer LT C400

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications

Technical

Gross volume:	383 litres
Temp. range:	-10°C to -45°C
Refrigerant:	R404a
Baskets:	1
Power consumption 24/h:	5,8 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage

Features

- Audible and visual alarm
- High/low alarm
- Battery backup (optional)
- Prepared for CO2 backup
- Castors
- Handle lock

Outer dimensions

Height:	885 mm
Width:	1560 mm
Depth:	695 mm
Gross weight:	95 kg
Net weight:	88 kg

Inner dimensions

Height:	635/385 mm
Width:	1190/210 mm
Depth:	440 mm

Product no.	N111004
-------------	---------



Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	24	264	N311220
			= 282	
75 mm	4	3	12	N311300
	8	24	192	N311320
			= 204	
100 mm	3	3	9	N311400
	5	24	120	N311420
			= 129	



Authorized Business Associate:

ANALYSIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811
Email: info@analysisindia.com • Web.: www.analysisindia.com