

Upright freezer LT U250

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications	
Technical	
Gross volume:	253 litres
Temp. range:	-10°C to -45°C
Refrigerant:	Nature R
Inner doors:	-
Adjustable shelves:	5
Power consumption 24/h:	5,5 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
Features	
- Temperature logger with USB port	
- Door alarm	
- High/low alarm	
- 48h battery backup (optional)	
- Castors	
Outer dimensions	
Height:	2055 mm
Width:	600 mm
Depth:	595 mm
Gross weight:	125 kg
Net weight:	112 kg
Inner dimensions	
Height:	1420 mm
Width:	420 mm
Depth:	370 mm
Product no.	N121002

Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	4 x 2 = 8	3	24	N321200
	5 x 2 = 10	12	120	N321210
			= 144	
75 mm	2 x 2 = 4	3	12	N321300
	3 x 2 = 6	12	72	N321310
			= 84	
100 mm	2 x 2 = 4	3	12	N321400
	2 x 2 = 4	12	48	N321400
			= 60	



Authorized Business Associate:

ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
 Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811
 Email: info@analysisindia.com • Web.: www.analysisindia.com