

# Chest freezer ULT C200

- Low cost of ownership
- Limited maintenance
- On-site service
- Easily installable & movable
- Eco-friendly, HCFC-free Nature R refrigerant
- Efficient cyclopentane insulation
- Low energy consumption
- ISO 9001, 14001 & 18001 OHSAS Certified

Specifications	
<b>Technical</b>	
Gross volume	198 litres
Internal storage capacity	189 litres
Temp. range	-60°C to -86°C
Temperature probes	1
Compressor	Embraco
Refrigerant	Nature R
Noise level	55 dB
Power consumption 24/h	6,0 kWh
Power consumption	2.190 kWh / year
Voltage	220V-240V
Frequency	50 Hz
<b>Features</b>	
Audible & visual alarm	YES
Castors with brakes	4
Handle lock	YES
Inner lids	2
Defrost	Manual
Battery Backup	50h
USB	NO
Data connection	NO
Digital controller with display	XR30CX
Data logging	NO
High/Low temperature alarm	YES
Analogue chart recorder	NO
Prepared for CO2 backup	YES
Prepared for external logging equipment	YES
<b>Dimensions</b>	
Height	850 mm
Height w. castors	890 mm
Width	920 mm
Depth	695 mm
Gross weight	69 kg
Net weight	64 kg
<b>Internal dimensions</b>	
Height	635/385 mm
Width	550/210 mm
Depth	440 mm
Insulation	80 mm
Insulation type	PUR
Internal humidity	<5 %
Product no.	N113002



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	9	99	N311220
			<b>= 117</b>	
75 mm	4	3	12	N311300
	8	9	72	N311320
			<b>= 84</b>	
100 mm	3	3	9	N311400
	5	9	45	N311420
			<b>= 54</b>	



Authorized Business Associate:

**ANALYSIS INSTRUMENTS LLP**

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)  
 Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811  
 Email: info@analysisindia.com • Web: www.analysisindia.com