

Chest freezer ULT C400

- Low cost of ownership
- Limited maintenance
- On-site service
- Easily installable & movable
- Eco-friendly, HCFC-free Nature R refrigerant
- Efficient cyclopentane insulation
- Low energy consumption
- ISO 9001, 14001 & 18001 OHSAS Certified

Specifications

Technical

Gross volume	383 litres
Internal storage capacity	368 litres
Temp. range	-60°C to -86°C
Temperature probes	1
Compressor	Embraco
Refrigerant	Nature R
Noise level	55 dB
Power consumption 24/h	8,10 kWh
Power consumption	2.957 kWh / year
Voltage	220V-240V
Frequency	50 Hz

Features

Audible & visual alarm	YES
Castors with brakes	4
Handle lock	YES
Inner lids	4
Defrost	Manual
Battery Backup	50h
USB	NO
Data connection	NO
Digital controller with display	XR30CX
Data logging	NO
High/Low temperature alarm	YES
Analogue chart recorder	NO
Prepared for CO2 backup	YES
Prepared for external logging equipment	YES

Dimensions

Height	850 mm
Height w. castors	890 mm
Width	1560 mm
Depth	695 mm
Gross weight	95 kg
Net weight	88 kg

Internal dimensions

Height	635/385 mm
Width	1190/210 mm
Depth	440 mm
Insulation	80 mm
Insulation type	PUR
Internal humidity	<5 %

Product no.	N113004
-------------	---------

ULT C400



Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	24	264	N311220
			= 282	
75 mm	4	3	12	N311300
	8	24	192	N311320
			= 204	
100 mm	3	3	9	N311400
	5	24	120	N311420
			= 129	



analysis®

JAS-ANZ



ISO 9001:2008
CERTIFIED

Authorized Business Associate:

ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
Mumbai - 400 002, INDIA • Tel.: +91-22-2205 2752 / 9820208811
Email: info@analysisindia.com • Web.: www.analysisindia.com