

# Chest freezer ULT C75



- Low cost of ownership
- Limited maintenance
- On-site service
- Easily installable & movable

- Eco-friendly, HCFC-free Nature R refrigerant
- Efficient cyclopentane insulation
- Low energy consumption
- ISO 9001, 14001 & 18001 OHSAS Certified

## Specifications

### Technical

Gross volume	74 litres
Internal storage capacity	71 litres
Temp. range	-60°C to -86°C
Temperature probes	1
Compressor	Embraco
Refrigerant	Nature R
Noise level	53 dB
Power consumption 24/h	5,10 kWh
Power consumption	1.862 kWh / year
Voltage	220V-240V
Frequency	50 Hz

### Features

Audible & visual alarm	YES
Castors with brakes	4
Handle lock	YES
Inner lids	1
Defrost	Manual
Battery Backup	50h
USB	NO
Data connection	NO
Digital controller with display	XR30CX
Data logging	NO
High/Low temperature alarm	YES
Analogue chart recorder	NO
Prepared for CO2 backup	YES
Prepared for external logging equipment	YES

### Dimensions

Height	895 mm
Height w. castors	935 mm
Width	550 mm
Depth	645 mm
Gross weight	57 kg
Net weight	53 kg

### Internal dimensions

Height	465 mm
Width	390 mm
Depth	390 mm
Insulation	80 mm
Insulation type	PUR
Internal humidity	<5 %

Product no.	N113001
-------------	---------

## ULT C75



## Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	8	5	40	N311210
75 mm	5	5	25	N311310
100 mm	4	5	20	N311410



Authorized Business Associate:

## ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)  
Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811  
Email: info@analisisindia.com • Web.: www.analisisindia.com