

Chest freezer XLT C150

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications

Technical

Gross volume:	140 litres
Temp. range:	-30°C to -60°C
Refrigerant:	Nature R
Baskets:	1
Power consumption 24/h:	3,5 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage

Features

- Audible and visual alarm
- High/low alarm
- Battery backup (optional)
- Prepared for CO2 backup
- Castors
- Handle lock

Outer dimensions

Height:	885 mm
Width:	720 mm
Depth:	695 mm
Gross weight:	59 kg
Net weight:	55 kg

Inner dimensions

Height:	635/385 mm
Width:	350/210 mm
Depth:	440 mm

Product no.	N112002
-------------	---------



Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	6	66	N311220
			= 84	
75 mm	4	3	12	N311300
	8	6	48	N311320
			= 60	
100 mm	3	3	9	N311400
	5	6	30	N311420
			= 39	



Authorized Business Associate:

ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811
Email: info@analysisindia.com • Web.: www.analysisindia.com