

Chest freezer XLT C300

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications	
Technical	
Gross volume:	296 litres
Temp. range:	-30°C to -60°C
Refrigerant:	Nature R
Baskets:	3
Power consumption 24/h:	6,2 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
Features	
- Audible and visual alarm	
- High/low alarm	
- Battery backup (optional)	
- Prepared for CO2 backup	
- Castors	
- Handle lock	
Outer dimensions	
Height:	885 mm
Width:	1260 mm
Depth:	695 mm
Gross weight:	82 kg
Net weight:	76 kg
Inner dimensions	
Height:	635/385 mm
Width:	890/210 mm
Depth:	440 mm
Product no.	N112003



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	18	198	N311220
			= 216	
75 mm	4	3	12	N311300
	8	18	144	N311320
			= 156	
100 mm	3	3	9	N311400
	5	18	90	N311420
			= 99	



Authorized Business Associate:

ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811
Email: info@analysisindia.com • Web.: www.analysisindia.com