

# Chest freezer XLT C400

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



## Specifications

### Technical

Gross volume:	383 litres
Temp. range:	-30°C to -60°C
Refrigerant:	Nature R
Baskets:	4
Power consumption 24/h:	7,7 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage

### Features

- Audible and visual alarm
- High/low alarm
- Battery backup (optional)
- Prepared for CO2 backup
- Castors
- Handle lock

### Outer dimensions

Height:	885 mm
Width:	1560 mm
Depth:	695 mm
Gross weight:	95 kg
Net weight:	88 kg

### Inner dimensions

Height:	635/385 mm
Width:	1190/210 mm
Depth:	440 mm

Product no.	N112004
-------------	---------



## Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	24	264	N311220
			<b>= 282</b>	
75 mm	4	3	12	N311300
	8	24	192	N311320
			<b>= 204</b>	
100 mm	3	3	9	N311400
	5	24	120	N311420
			<b>= 129</b>	



Authorized Business Associate:

## ANALYSIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)  
Mumbai - 400 002. INDIA • Tel.: +91-22-2205 2752 / 9820208811  
Email: info@analysisindia.com • Web.: www.analysisindia.com