

# Chest freezer XLT C500

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications	
<b>Technical</b>	
Gross volume:	495 litres
Temp. range:	-30°C to -60°C
Refrigerant:	Nature R
Baskets:	3
Power consumption 24/h:	6,0 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage
<b>Features</b>	
- Audible and visual alarm	
- High/low alarm	
- Battery backup (optional)	
- Prepared for CO2 backup	
- Legs	
- Handle lock	
<b>Outer dimensions</b>	
Height:	890 mm
Width:	1660 mm
Depth:	758 mm
Gross weight:	122 kg
Net weight:	107 kg
<b>Inner dimensions</b>	
Height:	640/390 mm
Width:	1290/210 mm
Depth:	500 mm
Product no.	N112005



Storage capacity				
Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	6	3	18	N311200
	11	27	297	N311220
			<b>= 315</b>	
75 mm	4	3	12	N311300
	8	27	216	N311320
			<b>= 228</b>	
100 mm	3	3	9	N311400
	5	27	135	N311420
			<b>= 144</b>	



Authorized Business Associate:

## ANALYSIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)  
 Mumbai - 400 002, INDIA • Tel.: +91-22-2205 2752 / 9820208811  
 Email: info@analysisindia.com • Web.: www.analysisindia.com