

Chest freezer XLT C75

- Low cost of ownership
- Limited maintenance
- On-site service
- Easy removable



Specifications

Technical

Gross volume:	74 litres
Temp. range:	-30°C to -60°C
Refrigerant:	Nature R
Baskets:	1
Power consumption 24/h:	4,0 kWh
Voltage:	220V-240V 50/60Hz
Compressor type:	Single stage

Features

- Audible and visual alarm
- High/low alarm
- Battery backup (optional)
- Prepared for CO2 backup
- Castors
- Handle lock

Outer dimensions

Height:	895 mm
Width:	550 mm
Depth:	645 mm
Gross weight:	57 kg
Net weight:	53 kg

Inner dimensions

Height:	465 mm
Width:	390 mm
Depth:	390 mm

Product no.	N112001
-------------	---------



Storage capacity

Box height	No. of boxes	No. of racks	Total no. of boxes	Product no.
50 mm	8	5	40	N311210
75 mm	5	5	25	N311310
100 mm	4	5	20	N311410



Authorized Business Associate:

ANALISIS INSTRUMENTS LLP

6, Sai Bhuvan, 3rd Floor, 141/145, S. G. Marg (Princess Street)
 Mumbai - 400 002, INDIA • Tel.: +91-22-2205 2752 / 9820208811
 Email: info@analisisindia.com • Web.: www.analisisindia.com